

pharmaxis	FM-361-04	Date Effective: 22-Nov-18	Page 1 of 1
BLISTER RELEASE			


Specification Number:	SP-623 - 08	Blister Lot Number:	M23-045
Description:	Bronchitol Blister Strip Australia and Europe.		

Checklist	Initial and Date
Production records have been reviewed and comply with requirements. (PIs, forms and associated records other than LBs)	JW 05 Oct 23
Room/Equipment Logbook reviewed, entries and cleaning comply with requirements.	LB-136 Primary Packaging Room – B1.022 - Logbook JW 05 Oct 23
Test results conform to specification	JW 05 Oct 23
Packaging materials released and within expiry at time of consumption	JW 05 Oct 23
Micro results have been entered in QC APQR Log	JW 05 Oct 23

		Record details	Initial and Date
Incident reports appropriately resolved (DR/NEPR/OOS/MTTR/CR/CAR/BMS OOL/BBM)		N/A.	JW 05 Oct 23
Retention sample number		223-208, 209, 210	JW 05 Oct 23
Capsule lots packed	0mg		
	5mg		
	10mg	N/A. 4 capsules 05 Oct 23	
	20mg		
	40mg	M23-039.	JW 05 Oct 23
Expiry date of capsules		23 Oct 23.	JW 05 Oct 23
Date of blister packing		8, 11 and 12-Sep-23	JW 05 Oct 23
Capsules have been released and were within expiry at time of consumption			JW 05 Oct 23
Batch number of blisters		LOT 23045	JW 05 Oct 23
Expiry date of blisters		EXP 07/2026	JW 05 Oct 23
No restrictions/conditions for release <ul style="list-style-type: none"> All tests required for all markets as per SP performed and comply No restrictions/conditions arising from incident reports No restrictions/conditions arising from restrictions/conditions on use of starting, packaging materials or intermediates. 		No restrictions	JW 05 Oct 23

Details of conditional release/Comments	
N/A. 4 capsules 05 Oct 23	
Quantity on paperwork (units)	27880
Quantity entered in Navision (units)	27880
Document quantity was derived from	P1-122-12, Section B.1.

Batch meets requirements for (tick box)	Full release	<input checked="" type="checkbox"/>	Conditional release	<input type="checkbox"/>
QA released by (sign and date):	J. Walsh 05 Oct 23			
Production record scanned, ELB-091 completed	J. Walsh 10 Oct 23			

Authorisation signature and date:	 05 Nov 18	Reference: SOP-295
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.		

Not entered on day of action JW 10 Oct 23

TB-002-03

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 1 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JK</i> 01 Sep 23		Lot Number: <i>M23-045</i>	

	Name	Position	Signature	Date
Author:	G.Velagala	Production Supervisor	<i>V. Velagala</i>	10 Dec 21
Checked by	Michael Lay	Production Operator	<i>Michael Lay</i>	10 Dec 21
Checked by (QA)	A. Dooley	QA Team Leader	<i>A. Dooley</i>	10 Dec 21
Authorised by:	S. Duffy	Production Manager	<i>S. Duffy</i>	10 Dec 21

1 PRODUCT DESCRIPTION - completed by Production Supervisor/Packaging Leader

1.1	Blister Description	Bronchitol Blister Strip Australia and Europe
	Specification Number (incl. version number)	SP-623-..... <i>09</i>
	Blister Lot Size (approx.)	* 40,000 30,000 Blisters <i>1215 SEP 23</i>
	Sample Plan (dependant on lot size).	Total blisters required (from spec): <u>155</u> Start (qty): <u>51</u> blisters Middle (qty): <u>53</u> blisters at each of following blister points: * 20,000 15,000 <i>1215 SEP 23</i> End (qty): <u>51</u> blisters
	Batch Number to be Debossed on Blister	<i>LOT 23045</i>
	Expiry Date to be Debossed on Blister	<i>EXP 07/2026</i>
	B1240 Control System recipe: (include version number)	4 Format 10 - <i>02</i>
	VisioChrom Control System recipe: (include version number)	VC_Bronchitol 10'S - <i>09</i>
	Filled Capsules Lot Number	<i>M23-039</i>
	Lid Foil Description & Specification Number	Bronchitol Foil Australia & Europe 20RR (SP-603)
	Clean required for the major equipment after batch (circle appropriate clean) as per SOP-271	<input checked="" type="radio"/> Major <input type="radio"/> Minor
COMPLETED BY/DATE: <i>JK</i> 05 SEP 23		

- **ENSURE SAFETY PRECAUTIONS ARE TAKEN AT ALL THE TIMES.**
- **FOLLOW THE STEPS IN ORDER AS SET OUT IN THE PI.**
- **REPORT ANY INCIDENTS TO THE SUPERVISOR IMMEDIATELY.**

2. The blister lot size is an estimate based on the output of the Ziegler machine process, the batch size is being changed to 30,000 and middle sample to 15,000 as the output from lot M23-031 only allows 30,000 blisters and 15,000

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 2 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JK 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

2 BILL OF MATERIALS

Material Name	Specification Number	Quantity per unit/Batch	Theoretical estimation of qty required	Lot Number	Expiry date
40mg Filled Printed Clear Capsules	SP-501	0.9g /Blister	27,000 g 36,000 *	<i>M23-034</i>	<i>23 OCT 23</i>
Blister Form Foil 160mm	SP-627	3g/Blister	90,000 g 120,000 *	<i>Tgn22-080</i> <i>N/A</i> <i>N/A</i>	<i>18 OCT 23</i>
Bronchitol Foil Australia and Europe – 20RR	SP-603	0.8 g/Blister	24,000 g 32,000 *	<i>Tgn22-047</i> <i>N/A</i> <i>N/A</i>	<i>10 MAY 25</i>
COMPLETED BY/DATE: <i>gg 07 Sep 23</i>					

*xx lot size had been adjusted. gg 19 Sep 23

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 3 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JR 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

3 ROOM & PRINCIPAL EQUIPMENT

3.1	<ul style="list-style-type: none"> Primary Packaging – 20RR (B1.022) Uhlmann B1240 Blister Thermoformer (Tag 00407) Calibrated Balances
-----	--

4 CLEANING CHECKS

4.1	Is this the first lot to be manufactured after a major clean? <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO If Yes, Check If more than 5 days have elapsed since the equipment was last sanitised, then: <ul style="list-style-type: none"> Re-sanitise all product contact surfaces with 70% IPA. Record re-sanitisation on the Equipment Cleaning Tag (FM-231) & in the primary packaging Log book (LB-136). Ensure all surfaces are dry. If No, this section is not applicable. Sign below when it is confirmed that the equipment is within the sanitisation expiry time. COMPLETED BY/DATE: <i>JR 07 Sep 23</i>
4.2	If it is the first lot since a major clean, then attach FM-297 to the processing instructions. Refer SOP-269. COMPLETED BY/DATE: <i>JR 07 Sep 23</i>
4.3	Attach FM-231 to the processing instructions. Refer SOP-269. COMPLETED BY/DATE: <i>JR 07 Sep 23</i>

5 ROOM CLEARANCE & EQUIPMENT SET-UP

5.1	Log in to the Primary Packaging Room (B1.022) - Logbook (LB-136). COMPLETED BY/DATE: <i>JR 07 Sep 23</i>
5.2	Ensure the room is clear of all unnecessary material – refer to SOP-142. Room clearance to be performed as per FM-391 by the second trained person (not involved in completing the above steps). COMPLETED BY/DATE: <i>JR 07 Sep 23</i> ROOM CLEARANCE CHECKED BY/DATE: <i>MQ 07 Sep 23</i>
5.3	Issue the foil. Check foils specification against Bill of Materials (Section 2). Ensure correct foils are issued and checked by a second trained person. Record foil details (directly from foil label – 1 label per roll of foil) and gross weight on FM-302. COMPLETED BY/DATE: <i>JR 07 Sep 23</i>
5.4	Set-up the Uhlmann B1240 Blister Thermoformer (Tag 00407) according to SOP-270. Record set-up on FM-301 (Uhlmann B1240 Format 10 Setup Checklist - 20RR). Whilst set-up is occurring, the process may continue.

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 4 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>PK 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

6 FILLED CAPSULES

6.1	<p>If the blister packaging lot is a continuation of a campaign using the same filled capsules as the previous lot, then record the previous blister packaging lot details:</p> <p>Previous Blister Lot Number: <i>N/A</i></p> <p>Previous Blister Description: <i>N/A</i></p> <p>COMPLETED BY/DATE: <i>gg 07 SEP 23</i></p>										
6.2	<p>Complete the following table:</p> <table border="1"> <tr> <td>Capsule Description</td> <td>40mg Filled Printed Clear Capsules</td> </tr> <tr> <td>SP Number</td> <td>SP-501</td> </tr> <tr> <td>Lot Number</td> <td><i>M23 -039</i></td> </tr> <tr> <td>Expiry Date (ensure capsules are within expiry)</td> <td><i>23 OCT 23</i></td> </tr> <tr> <td>Bag Numbers Issued</td> <td><i>1 - 10</i></td> </tr> </table> <p>COMPLETED BY/DATE: <i>gg 07 SEP 23</i></p>	Capsule Description	40mg Filled Printed Clear Capsules	SP Number	SP-501	Lot Number	<i>M23 -039</i>	Expiry Date (ensure capsules are within expiry)	<i>23 OCT 23</i>	Bag Numbers Issued	<i>1 - 10</i>
Capsule Description	40mg Filled Printed Clear Capsules										
SP Number	SP-501										
Lot Number	<i>M23 -039</i>										
Expiry Date (ensure capsules are within expiry)	<i>23 OCT 23</i>										
Bag Numbers Issued	<i>1 - 10</i>										
6.3	<p>Label the following sterile plastic bags. Refer SOP-037:</p> <p>Reject Capsules (Label "Reject" & with filled capsule lot details). Record the tare weight on label.</p> <p>Waste Foil (Label "Reject" with blister lot details).</p> <p>Filled Blisters to be Stripped (Label "WIP" with blister lot details).</p> <p>COMPLETED BY/DATE: <i>gg 08 Sep 23</i></p>										
6.4	<p>Set-up of Uhlmann Blisterpacker is complete once the blister integrity test is complete and a 'Pass' result was achieved & recorded.</p> <ul style="list-style-type: none"> Attach the completed Format Setup Checklist (FM-301) to the PI. <p>COMPLETED BY/DATE: <i>gg 08 SEP 23</i></p>										

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 5 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JK 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

7 RECIPE & MACHINE DETAILS

7.1	<p>Switch on the B1240 as described in SOP-269 & login to the B1240 Control System & VisioChrom Control System.</p> <p><u>Record recipes selected and ensure in compliance with Section 1.1</u></p> <p>B1240 Control System recipe (incl. version number): 4 FORMAT 10-02</p> <p>VisioChrom Control System recipe (incl. version number): VC-Bronchitol 105-09</p> <p>Allow sealing plate temperature to increase and settle within range. If it does not settle within range, contact Supervisor/Packaging Leader.</p> <p>Record the sealing plate temperature & machine speed:</p> <table border="1" style="width: 100%;"> <tr> <td style="width: 40%;">Sealing Plate Temperature:</td> <td style="width: 20%; text-align: center;">179 °C</td> <td style="width: 40%; text-align: center;">179-181°C</td> </tr> <tr> <td>Machine Speed:</td> <td style="text-align: center;">40 Cycles/min</td> <td style="text-align: center;">SET</td> </tr> </table> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>	Sealing Plate Temperature:	179 °C	179-181°C	Machine Speed:	40 Cycles/min	SET
Sealing Plate Temperature:	179 °C	179-181°C					
Machine Speed:	40 Cycles/min	SET					

8 ROOM CONDITIONS

8.1	<p>Check for EMS alarms. If no alarms are occurring on the EMS regarding the pressure differentials, temperatures and relative humidity in Room B1.022, then it is accepted that all values are within their acceptable ranges for production to occur.</p> <p>If an alarm is present, do not proceed with the batch until the conditions are restored to within acceptable limits (no alarms) or a Deviation report (FM-016) has been raised and approved.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>
-----	--

9 PREPARATION FOR START UP

9.1	<ul style="list-style-type: none"> Re-position the Compensation Pendulum to the running position. Place the "Waste Foil" bag under the machine reject outlet. <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>										
9.2	<p>Start the machine and perform batch challenge testing for 'Reject Mechanism' as detailed in SOP-269 and record the results below when performed:</p> <table border="1" style="width: 100%;"> <thead> <tr> <th>Test</th> <th>Challenge Test</th> <th>Expected Results</th> <th>Lane 1 (Passed=Yes/No)</th> <th>Lane 2 (Passed=Yes/No)</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Blister with Empty pockets</td> <td>Reject</td> <td style="text-align: center;">yes</td> <td style="text-align: center;">yes</td> </tr> </tbody> </table> <p>Do not proceed if the test does not comply, report minor incident and inform production supervisor immediately. Sign/date below when it is confirmed that test(s) complies.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>	Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)	1	Blister with Empty pockets	Reject	yes	yes
Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)							
1	Blister with Empty pockets	Reject	yes	yes							

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 6 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date: <i>JK: 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

9.3	<p>Label each hopper (in pouches provided) with a WIP label showing filled capsule details. Pour the capsules directly from the bag(s) into the hoppers.</p> <p>A small amount of capsules should be poured from the filled capsule bag into a sanitised 'Hand Filling Capsules' container. Label the container with a WIP label referencing 'Hand Filling Capsules' & filled capsule details.</p> <p>Record the addition time/date, balance details, capsule details & bag weights on FM-302 as capsules are introduced into the hoppers. All weights must be checked by a 2nd operator.</p> <table border="1"> <thead> <tr> <th>Filled Capsule Description</th> <th>Fill Pockets By</th> </tr> </thead> <tbody> <tr> <td>40mg Filled Printed Clear Capsules</td> <td>SIMTAP Feeder Hoppers</td> </tr> </tbody> </table> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>	Filled Capsule Description	Fill Pockets By	40mg Filled Printed Clear Capsules	SIMTAP Feeder Hoppers											
Filled Capsule Description	Fill Pockets By															
40mg Filled Printed Clear Capsules	SIMTAP Feeder Hoppers															
9.4	<p>Train the Visio Chrom camera as per SOP-269.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>															
9.5	<ul style="list-style-type: none"> Set the SIMTAP Feeder switch to position 1 (manual) & fill the Sorting Plate until capsules stop feeding. Switch the SIMTAP Feeder switch to position 2 (automatic). Press reset & then start the machine. <p><i>Ref: CM-23-003- 08 Sep 23</i></p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>															
9.6	<p>Perform start of the batch challenge testing for Visio Chrom Camera as detailed in SOP-269 and record the results below when performed:</p> <table border="1"> <thead> <tr> <th>Test</th> <th>Challenge Test</th> <th>Expected Results</th> <th>Lane 1 (Passed=Yes/No)</th> <th>Lane 2 (Passed=Yes/No)</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Empty Pocket</td> <td>Reject</td> <td><i>Yes</i></td> <td><i>Yes</i></td> </tr> <tr> <td>2</td> <td>Good Capsules in all pockets</td> <td>Accept</td> <td><i>Yes</i></td> <td><i>Yes</i></td> </tr> </tbody> </table> <p>Do not proceed if the test (s) does not comply record minor incident and inform production supervisor immediately. Sign/date below when it is confirmed that test(s) complies.</p> <p>Place capsules used for camera verification and blisters produced in the corresponding reject bag.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>	Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)	1	Empty Pocket	Reject	<i>Yes</i>	<i>Yes</i>	2	Good Capsules in all pockets	Accept	<i>Yes</i>	<i>Yes</i>
Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)												
1	Empty Pocket	Reject	<i>Yes</i>	<i>Yes</i>												
2	Good Capsules in all pockets	Accept	<i>Yes</i>	<i>Yes</i>												
9.7	<p>Run the machine until filled blisters reach the punching station and stop immediately.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>															
9.8	<p>Reset statistical data as detailed in SOP-269.</p> <p>COMPLETED BY/DATE: <i>JK 08 Sep 23</i></p>															
9.9	<p>Line clearance must be performed as per SOP-142 and recorded on FM-304 (Line Clearance Checklist for Blister Packaging Processes at Pharmaxis).</p> <p>LINE CLEARANCE PERFORMED BY/DATE: <i>AL 08 Sep 23</i></p>															

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 7 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JR</i> <i>01 Sep 23</i>		Lot Number: <i>M23-045</i>	

10 BATCH MANUFACTURE

10.1	Commence batch blister packing. Record training of the camera, in-process checks, leak test sampling, weights of raw / WIP materials & machine adjustments in FM-302. Refer to the table below & SOP-269.																																													
	<table border="1"> <thead> <tr> <th>Check-point</th> <th>Train Camera</th> <th>In-Process Check</th> <th>Leak Test (6 blisters – 3 from each lane)</th> <th>Sealing Plate temperature</th> </tr> </thead> <tbody> <tr> <td>Start of Production / Day</td> <td>Yes</td> <td>Yes</td> <td>Yes</td> <td>Yes</td> </tr> <tr> <td>Hourly Check</td> <td>No</td> <td>Yes</td> <td>Yes</td> <td>No</td> </tr> <tr> <td>Lid Foil - end of roll</td> <td>No</td> <td>Yes</td> <td>Yes</td> <td>No</td> </tr> <tr> <td>Lid Foil - start of roll</td> <td>No</td> <td>Yes</td> <td>Yes</td> <td>No</td> </tr> <tr> <td>Form Foil Change</td> <td>No</td> <td>Yes</td> <td>Yes</td> <td>No</td> </tr> <tr> <td>After break</td> <td>No</td> <td>Yes</td> <td>No</td> <td>No</td> </tr> <tr> <td>End of Production / Day</td> <td>No</td> <td>Yes</td> <td>Yes</td> <td>No</td> </tr> <tr> <td>QC samples</td> <td colspan="4">As described in section 1. Samples must be taken equally from each lane. Label and handle QC samples as described in SOP-269.</td> </tr> </tbody> </table>	Check-point	Train Camera	In-Process Check	Leak Test (6 blisters – 3 from each lane)	Sealing Plate temperature	Start of Production / Day	Yes	Yes	Yes	Yes	Hourly Check	No	Yes	Yes	No	Lid Foil - end of roll	No	Yes	Yes	No	Lid Foil - start of roll	No	Yes	Yes	No	Form Foil Change	No	Yes	Yes	No	After break	No	Yes	No	No	End of Production / Day	No	Yes	Yes	No	QC samples	As described in section 1. Samples must be taken equally from each lane. Label and handle QC samples as described in SOP-269.			
Check-point	Train Camera	In-Process Check	Leak Test (6 blisters – 3 from each lane)	Sealing Plate temperature																																										
Start of Production / Day	Yes	Yes	Yes	Yes																																										
Hourly Check	No	Yes	Yes	No																																										
Lid Foil - end of roll	No	Yes	Yes	No																																										
Lid Foil - start of roll	No	Yes	Yes	No																																										
Form Foil Change	No	Yes	Yes	No																																										
After break	No	Yes	No	No																																										
End of Production / Day	No	Yes	Yes	No																																										
QC samples	As described in section 1. Samples must be taken equally from each lane. Label and handle QC samples as described in SOP-269.																																													
	COMPLETED BY/DATE: <i>JR 12 SEP 23</i>																																													
10.2	When the requisite number of blisters has been prepared and the amount obtained in secondary is ensured: Switch the SIMTAP feeder to position zero <i>Ref. CM-23-003-D1 08/12 Sep 23</i> Run off the good blisters to secondary until the SIMTAP feeder is empty Stop the machine as per SOP-269 COMPLETED BY/DATE: <i>JR 12 SEP 23</i>																																													

11 ACTIONS AT COMPLETION OF LOT BLISTER PACKING

11.1	Start the machine and perform end of the batch challenge test for 'Reject Mechanism' as detailed in SOP-269 and record the results below when performed:										
	<table border="1"> <thead> <tr> <th>Test</th> <th>Challenge Test</th> <th>Expected Results</th> <th>Lane 1 (Passed=Yes/No)</th> <th>Lane 2 (Passed=Yes/No)</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Blister with Empty pockets</td> <td>Reject</td> <td><i>yes</i></td> <td><i>yes</i></td> </tr> </tbody> </table>	Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)	1	Blister with Empty pockets	Reject	<i>yes</i>	<i>yes</i>
Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)							
1	Blister with Empty pockets	Reject	<i>yes</i>	<i>yes</i>							
	Do not proceed if the test does not comply record minor incident and inform production supervisor immediately. Sign/date below when it is confirmed that test(s) complies. COMPLETED BY/DATE: <i>JR 12 SEP 23</i>										
11.2	Perform end of the batch challenge testing for Visio Chrom Camera as detailed in SOP-269 and record the results below when performed:										

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 8 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JP 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

Test	Challenge Test	Expected Results	Lane 1 (Passed=Yes/No)	Lane 2 (Passed=Yes/No)
1	Empty Pocket	Reject	<i>yes</i>	<i>yes</i>
2	Good Capsules in all pockets	Accept	<i>yes</i>	<i>yes</i>

Do not proceed if the test (s) does not comply record minor incident and inform production supervisor immediately. Sign/date below when it is confirmed that test(s) complies.

Place capsules used for camera verification and blisters produced in the corresponding reject bag.

COMPLETED BY/DATE: *JP 12 Sep 23*

11.3 Place all remaining capsules in the SIMTAP Distributor Plate, capsules remaining in the 'Hand Filling' container and blisters produced from end of batch challenge testing in the corresponding "Reject" labelled bag, if applicable.

COMPLETED BY/DATE: *JP 12 Sep 23*

11.4 Save the lot Statistical Data to a removable storage device.
The '.pdf' file should have the following format: LOT NUMBER – DATE.
Log out of the B1240 and VisioChrom control system – refer to SOP-269.

COMPLETED BY/DATE: *JP 12 Sep 23*

11.5 Remove unused lid foil & form foil from the machine. Unused foil is to be weighed & checked by trained operators. Record tare weights on FM-302.
Return unused foil to the Warehouse.

COMPLETED BY/DATE: *JP 12 Sep 23*

11.6 At the end of the batch, unused capsules (if any) remaining in the hoppers may be recovered. Label each recovered capsule bag with a WIP label, weigh empty and record tare weight & the description on WIP label as "Recovered 40mg Filled Printed Clear Capsules – 20RR". Weigh the recovered (i.e. Unused) capsule bag and record the Gross Weight below. Calculate the Net Weight.
Weights and calculations must be checked by a 2nd Operator. *JP 12 Sep 23*

Balance Used (Tag): _____

Description	Gross Weight	Tare Weight	Net Weight
Recovered Unused Capsules	g	g	g

COMPLETED BY/DATE: _____

CHECKED BY/DATE: _____

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 9 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <u>JK 01 Sep 23</u>		Lot Number: <u>M23-045</u>	

11.7	Strip the capsules from the "Filled Blisters to be Stripped" blisters and place in the 'Reject Capsules' bag. COMPLETED BY/DATE: <u>gg 13 SEP 23</u>								
11.8	Weigh the Reject Capsule bag. All weights must be checked by a 2nd Operator. Balance Used (Tag): <u>00011</u> <table border="1"> <thead> <tr> <th>Description</th> <th>Gross Weight</th> <th>Tare Weight</th> <th>Net Weight</th> </tr> </thead> <tbody> <tr> <td>Reject Capsules</td> <td>2304.7 g</td> <td>42.3 g</td> <td>2262.4 g</td> </tr> </tbody> </table> COMPLETED BY/DATE: <u>gg 13 Sep 23</u> CHECKED BY/DATE: <u>gg 13 Sep 23</u>	Description	Gross Weight	Tare Weight	Net Weight	Reject Capsules	2304.7 g	42.3 g	2262.4 g
Description	Gross Weight	Tare Weight	Net Weight						
Reject Capsules	2304.7 g	42.3 g	2262.4 g						
11.9	Enter consumption of materials into all associated inventory forms. Tick (✓) when the values for each material type have been entered. <table border="1"> <tr> <td>Capsules</td> <td><input checked="" type="checkbox"/></td> </tr> <tr> <td>Lid Foil</td> <td><input checked="" type="checkbox"/></td> </tr> <tr> <td>Blister Form Foil</td> <td><input checked="" type="checkbox"/></td> </tr> </table> Also ensure the correct values have been entered on all forms. COMPLETED BY/DATE: <u>gg 13 SEP 23</u> CHECKED BY/DATE: <u>gg 13 Sep 23</u>	Capsules	<input checked="" type="checkbox"/>	Lid Foil	<input checked="" type="checkbox"/>	Blister Form Foil	<input checked="" type="checkbox"/>		
Capsules	<input checked="" type="checkbox"/>								
Lid Foil	<input checked="" type="checkbox"/>								
Blister Form Foil	<input checked="" type="checkbox"/>								

12 MACHINE SHUTDOWN & CLEANING

12.1	Logout of the VisioChrom and Blister Thermoformer HMIs as described in SOP-269. Turn off main switch at the rear of the machine. The air supply automatically cuts off. Perform cleaning described in Sect 1.1 of PI. Type of cleaning performed: <u>Major</u> Note: Ensure that a major clean is done for the used format parts if the batch is the final lot of the campaign. Otherwise minor clean is to be completed. Contact Production Supervisor to confirm, as required. COMPLETED BY/DATE: <u>gg 14 Sep 23</u>
------	---

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 10 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date <i>JK 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

13 CONSUMPTION & RECONCILIATION

13.1	Calculate the total consumption of the following materials from FM-302. <table border="1"> <thead> <tr> <th>Material</th> <th>Spec. Number</th> <th>Lot Number(s)</th> <th>Total Quantity Used per Lot</th> </tr> </thead> <tbody> <tr> <td>40mg Filled Printed Clear Capsules</td> <td>SP-501</td> <td><i>M23-039</i></td> <td> Total issued on FM-302 – Recovered Capsules <i>gg 105023</i> 27693 + <i>27043.1</i> g </td> </tr> <tr> <td rowspan="3">Blister Form Foil 160mm (Record different lot numbers separately)</td> <td rowspan="3">SP-627</td> <td><i>Ign22-080</i></td> <td><i>77753.2</i> g</td> </tr> <tr> <td><i>N/A</i></td> <td><i>N/A</i> g</td> </tr> <tr> <td><i>N/A</i></td> <td><i>N/A</i> g</td> </tr> <tr> <td rowspan="3">Bronchitol Foil Australia and Europe – 20RR (Record different lot numbers separately)</td> <td rowspan="3">SP-603</td> <td><i>Ign22-041</i></td> <td><i>23206.5</i> g</td> </tr> <tr> <td><i>N/A</i></td> <td><i>N/A</i> g</td> </tr> <tr> <td><i>N/A</i></td> <td><i>N/A</i> g</td> </tr> </tbody> </table>	Material	Spec. Number	Lot Number(s)	Total Quantity Used per Lot	40mg Filled Printed Clear Capsules	SP-501	<i>M23-039</i>	Total issued on FM-302 – Recovered Capsules <i>gg 105023</i> 27693 + <i>27043.1</i> g	Blister Form Foil 160mm (Record different lot numbers separately)	SP-627	<i>Ign22-080</i>	<i>77753.2</i> g	<i>N/A</i>	<i>N/A</i> g	<i>N/A</i>	<i>N/A</i> g	Bronchitol Foil Australia and Europe – 20RR (Record different lot numbers separately)	SP-603	<i>Ign22-041</i>	<i>23206.5</i> g	<i>N/A</i>	<i>N/A</i> g	<i>N/A</i>	<i>N/A</i> g
Material	Spec. Number	Lot Number(s)	Total Quantity Used per Lot																						
40mg Filled Printed Clear Capsules	SP-501	<i>M23-039</i>	Total issued on FM-302 – Recovered Capsules <i>gg 105023</i> 27693 + <i>27043.1</i> g																						
Blister Form Foil 160mm (Record different lot numbers separately)	SP-627	<i>Ign22-080</i>	<i>77753.2</i> g																						
		<i>N/A</i>	<i>N/A</i> g																						
		<i>N/A</i>	<i>N/A</i> g																						
Bronchitol Foil Australia and Europe – 20RR (Record different lot numbers separately)	SP-603	<i>Ign22-041</i>	<i>23206.5</i> g																						
		<i>N/A</i>	<i>N/A</i> g																						
		<i>N/A</i>	<i>N/A</i> g																						
COMPLETED BY/DATE: <i>gg 13 SEP 23</i>																									
CHECKED BY/DATE: <i>gg 13 Sep 23</i>																									
13.2	Print a copy of the Statistical data file. Write Lot number and sign/date the printout. Attach to PI. COMPLETED BY/DATE: <i>gg 13 Sep 23</i>																								
13.3	At the end of each lot, complete the capsule reconciliation for the capsules used on MES-034. If reconciliation limit is exceeded, initiate a Deviation Report (as per SOP-030). Also record output (Blisters Transferred to Secondary Packaging). <table border="1"> <tr> <td>Output:</td> <td>Quantity Blisters (EA)</td> </tr> <tr> <td></td> <td><i>27880</i></td> </tr> </table>	Output:	Quantity Blisters (EA)		<i>27880</i>																				
Output:	Quantity Blisters (EA)																								
	<i>27880</i>																								
COMPLETED BY/DATE: <i>gg 13 Sep 23</i>																									

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 11 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date	J 01 Sep 23	Lot Number:	M23-045

14 MINOR INCIDENT REPORT

Record any observation made during the batch manufacture – refer to SOP-I142 for minor incident definition, identification and subsequent action.

Alert the Production Supervisor immediately. Action!

Date/Time	Operator Initials	Incident/Observation – brief description/actions	Supervisor Evaluation – root cause, product impact, immediate action	Product Quality Impact – If N – provide brief justification If Y – reference deviation report.	Supervisor Approval – Sign/Date
08 Sep 23	JB	Start of the batch still the same issued of the VC-concerning any defect that have one and some with Jemazation over to produce and increase	Both issues are known to production and engineering the current reject rate is same as previous runs which is acceptable for the current run.	No quality impact on final product: only performance issue and production all defects will be detected and reject during the run.	13 SEP 23
N/A	N/A	regard vial: Eng, many pri used informed.	N/A	N/A	N/A
				N/A 9/13 Sep 23	

pharmaxis	PI-122-12	Date Effective: 23-Dec-21	Page 12 of 12
PROCESSING INSTRUCTIONS: UHLMANN B1240 BRONCHITOL BLISTER STRIP AU/EU (SP-623)			
Issued By/Date: <i>[Signature]</i> 01 Sep 23		Lot Number: M23-045	

15 NAVISION & TRENDING UPDATE

15.1	Perform Navision Production Journal transactions as per SOP-165 and enter Trending as per SOP-241 COMPLETED BY/DATE: <i>[Signature]</i> 14 Sep 23
------	---

16 BATCH REVIEW

16.1	Verify Navision Production Journal transactions as per SOP-165. COMPLETED BY/DATE: <i>[Signature]</i> 15 SEP 23		
16.2	Ensure production journal is posted & trending is completed accurately. COMPLETED BY/DATE: <i>[Signature]</i> 19 Sep 23		
16.3	Record any Associated Deviation, CR or CAR – to be completed by Production staff <table border="1"> <tr> <td>Reports</td> <td>N/A</td> </tr> </table> COMPLETED BY/DATE: <i>[Signature]</i> 22 Sep 23	Reports	N/A
Reports	N/A		
16.4	The Process Instructions & the associated documents have been reviewed as detailed in SOP-291 COMPLETED BY/DATE: <i>[Signature]</i> 22 Sep 23		

17 APPROVAL BY PRODUCTION MANAGER/PRODUCTION SUPERVISOR

These Processing Instructions & associated documents have been reviewed & the lot produced in compliance with GMP.

Print Name: *Larry Reed* Sign/Date: *[Signature]* 05 Oct 23

18 ASSOCIATED DOCUMENTS AND QA REVIEW


Include details of any associated documents

Associated Document	Document Number	Initial/Date
EMS OOL	N/A	<i>[Signature]</i> 05 Oct 23
CAR Number(s)	N/A	<i>[Signature]</i> 05 Oct 23
DR Number(s)	N/A	<i>[Signature]</i> 05 Oct 23
NCPR Number(s)	N/A	<i>[Signature]</i> 05 Oct 23
CR Number(s)	N/A	<i>[Signature]</i> 05 Oct 23
OOS Number(s)	N/A	<i>[Signature]</i> 05 Oct 23

BATCH DOCUMENT QA REVIEWED BY

Print Name: *Subin Wat* Sign/Date: *[Signature]* 05 Oct 23


Lot Number: M23-045

Authorisation signature and date:  16 Mar 21	Reference: SOP-269	62X J912SEP23
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.		
TF-015-03		

Lot Number:	M23-045
-------------	---------

TE-015-03

Lot Number: M23-045

Authorisation signature and date:  16 Mar 21	Reference: SOP-269
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-015-03	

Lot Number: M23-045

94 OK Sep 23

Capsule Type (mg): _____; Capsule Lot No. _____; Expiry Date: _____

Capsule Type (mg): _____; **Capsule Lot No.** _____; **Expiry Date:** _____

TE-015-03

pharmaxis

FM-302-11

Date Effective:
30-Mar-21

Page 5 of 6

PROCESS CONTROL RECORD FOR BLISTER-PACKING PROCESS AT PHARMAXIS

Issued By/Date: *JK* 01 Sep 23

Lot Number: M23-045

LID FOIL - record details and weights used in table below.

Balance Tag #: 00011

Completed by/date *JK* 07 Sep 23

Date	Lot No.	Artwork Number/Version	Expiry Date	Container No. __ of __ (foil)	Gross Weight (g)	Completed by (Initials)	Checked by (Initials)	Date/Time (Immediately prior to load onto the machine)	Tare Weight (g)	Net Weight (g)	Completed by (Initials)	Checked by (Initials)
07 SEP 23	Ign22-041	SP603/02	19 MAY 25	10 OF 13	12410.	<i>JK</i>	<i>JK</i>	07 Sep 23 10:53	178.2	12231.8	<i>JK</i>	<i>JK</i>
07 Sep 23	Ign22-041	SP602/02	19 MAY 25	4 OF 13	13285.	<i>JK</i>	<i>JK</i>	11 Sep 23 10:07	2310.3	10974.7	<i>JK</i>	<i>JK</i>
<i>N/A</i> <i>JK</i> 13 SEP 23												
									Total (g)	23206.5		

pharmaxis

FM-302-11

Date Effective:
30-Mar-21

Page 6 of 6

PROCESS CONTROL RECORD FOR BLISTER-PACKING PROCESS AT PHARMAXIS

Issued By/Date: *JK* 01 Sep 23

Lot Number: M23-045

FORM FOIL - record details and weights used in table below.

Balance Tag #: 00011

Completed by/date: *gg* 07 Sep 23

Date	Lot No.	Artwork Number/Version	Expiry Date	Container No. __ of __ (foil)	Gross Weight (g)	Completed by (Initials)	Checked by (Initials)	Date/Time (Immediately prior to load onto the machine)	Tare Weight (g)	Net Weight (g)	Completed by (Initials)	Checked by (Initials)
07 Sep 23	Ign22-080	N/A	08 OCT 23	12 4 OF 28 <i>28 07 Sep 23</i>	4536.1	<i>gg</i>	<i>[initials]</i>	07 Sep 23 10:56	516.0	4020.1	<i>gg</i>	<i>[initials]</i>
07 Sep 23	Ign22-080	N/A	08 OCT 23	4 OF 28	28407.	<i>gg</i>	<i>[initials]</i>	08 Sep 23 10:24	516.0	27891.	<i>gg</i>	<i>[initials]</i>
08 Sep 23	Ign22-080	N/A	08 OCT 23	15 OF 28	28258.	<i>gg</i>	<i>[initials]</i>	11 Sep 23 09:03	515.9	27742.1	<i>gg</i>	<i>[initials]</i>
11 Sep 23	Ign22-080	N/A	08 OCT 23	1 OF 28	28471.	<i>gg</i>	<i>[initials]</i>	11 Sep 23 12:51	10371.	18100.	<i>gg</i>	<i>[initials]</i>
<div><div></div><div>N/A</div><div>gg 13 Sep 23</div></div>												
Total (g)										77753.2		

MES-034-05**Reconciliation of Filled Capsules After Blisterpacking on the Uhlmann B1240**

Modified 10-Jul-15

Excel 2016

Reference: SOP-269

Capsule Lot Number:	M23-039
Capsule Description:	40mg Filled printed Clear Capsules
Capsule Specification Number:	SP-501
Average Gross Capsule Weight (g):	0.08815 (from Encapsulation MES-031)

Blister Lot Number:	M23-045
Number of Capsules per Blister:	10 (from PI)
Weight of Capsules Issued (g):	27,093.1 (from FM-302)

	Blisters (Count)	Capsule Weight (g)	
Capsules Issued:		27,093.1	
Reject Capsules:		2,262.4	(from PI)
Recovered Capsules:		0.0	(from PI)
Leak Test Sample Blisters:	96	84.6	(only leak test sample blisters which contain capsules, from FM-302)
QA Sample Blisters:	158	139.3	(from FM-302)
Good Blisters:	28,134	24,800.1	(from Uhlmann .pdf)
Total Capsules Accounted for:		27,062.5	

Blister Lot Number:	Fill-in
Number of Capsules per Blister:	Fill-in (from PI)
Weight of Capsules Issued (g):	Fill-in (from FM-302)

	Blisters (Count)	Capsule Weight (g)	
Capsules Issued:		0.0	
Reject Capsules:		Fill-in	(from PI)
Recovered Capsules:		Fill-in	(from PI)
Leak Test Sample Blisters:	Fill-in	0.0	(only leak test sample blisters which contain capsules, from FM-302)
QA Sample Blisters:	Fill-in	0.0	(from FM-302)
Good Blisters:	Fill-in	0.0	(from Uhlmann .pdf)
Total Capsules Accounted for:		0.0	

Blister Lot Number:	Fill-in
Number of Capsules per Blister:	Fill-in (from PI)
Weight of Capsules Issued (g):	Fill-in (from FM-302)

	Blisters (Count)	Capsule Weight (g)	
Capsules Issued:		0.0	
Reject Capsules:		Fill-in	(from PI)
Recovered Capsules:		Fill-in	(from PI)
Leak Test Sample Blisters:	Fill-in	0.0	(only leak test sample blisters which contain capsules, from FM-302)
QA Sample Blisters:	Fill-in	0.0	(from FM-302)
Good Blisters:	Fill-in	0.0	(from Uhlmann .pdf)
Total Capsules Accounted for:		0.0	

Blister Lot Number:	Fill-in
Number of Capsules per Blister:	Fill-in (from PI)
Weight of Capsules Issued (g):	Fill-in (from FM-302)

	Blisters (Count)	Capsule Weight (g)	
Capsules Issued:		0.0	
Reject Capsules:		Fill-in	(from PI)
Recovered Capsules:		Fill-in	(from PI)
Leak Test Sample Blisters:	Fill-in	0.0	(only leak test sample blisters which contain capsules, from FM-302)
QA Sample Blisters:	Fill-in	0.0	(from FM-302)
Good Blisters:	Fill-in	0.0	(from Uhlmann .pdf)
Total Capsules Accounted for:		0.0	

RECONCILIATION OF CAPSULES - UHLMANN

Total Weight of Capsules Issued (g):	27,093.1
Total Weight of Capsules Accounted for (g):	27,062.5
% Reconciliation:	99.9 (Limits: 98.0-102.0%)
Result:	Pass (Pass / Fail)

BLISTERS TO SECONDARY PACKAGING ROOM

Blister Lot Number:	M23-045	N/A	N/A	N/A
Blisters Transferred to Secondary Packaging Room:	27,880	N/A	N/A	N/A

Completed by/date: gg 13 Sep 23Checked by/date: AD 13 Sep 23

Refer to any DR/CAR/CR/NCPR etc.

ROOM CLEARANCE CHECKLIST

Lot Number: M23-045

Issued By/Date: [Signature] 01 Sep 23

ROOM CLEARANCE CHECKS

	✓ each step once complete
Correct area, is in use (as specified in the process instructions).	✓
The outer door of the room is status labelled.	✓
All the documentation, products and materials from the previous production operation have been removed from the room.	✓
Ensure the following: <ul style="list-style-type: none">• Equipment is appropriate to use - cleaning completed as required. Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>• Equipment is clean & dry (visible inspection only). Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>• Sanitisation is within the time frame as specified in the PI. Yes <input type="checkbox"/> N/A <input checked="" type="checkbox"/>	✓
Cleaning tags are attached to the PI and ensure that the major clean has been approved (if applicable).	✓
All the steps up to the room clearance section are completed in the corresponding PI.	✓
All the logbooks involved are checked for completeness (e.g., Login, batch details, minor clean/daily clean performed).	✓
The Balances involved are calibrated and the daily check(s) performed.	✓
No alarms on EMS system (ensure alarm and strobe lighting not initiated).	✓
Only the Bill of Materials listed in the PI are present in the room and are within expiry date.	✓
Ensure that the details recorded in PI for the bill of materials are accurate.	✓
Comments: B1022	
If any of the checks above does not comply, do not proceed with Room Clearance.	
ROOM CLEARANCE COMPLETED BY/DATE: [Signature] 07 Sep 23	

Sign & Date in the PI & in the corresponding logbook upon successful completion of the Room Clearance.

Authorisation signature and date:

[Signature] 08 Jun 22

Reference: SOP-142

Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.

TE-015-03

pharmaxis

FM-304-08

Date Effective:
10-May-2023

Page 1 of 1

LINE CLEARANCE CHECKLIST FOR BLISTER PACKAGING PROCESSES AT PHARMAXIS

Lot Number: M23-045

Issued By/Date:

J 01 Sep 23

Product – Blister Description	Lot Number	Specification Number and Version
<i>Band-aid Blister Strip Australia and Europe</i>	<i>M23-045</i>	<i>SP-623-08</i>
Line Clearance Checks		✓ each step once complete
All processing instructions, associated forms, checklist etc. are present, issued by QA and correct Lot Number identified on each.		
No traces of documentation/materials from previous batch in the area.		
No rubbish visible in surrounding area. No tools present from set-up.		
No alarms on EMS system (ensure alarm and strobe lighting not initiated).		
Room Logbook complete (LB-136).		
Equipment is clean and dry (visual inspection only).		
Cleaning records (FM-231 and FM-297) complete, removed from machine & attached to PI.		
Daily balance check(s) performed and recorded in associated logbook(s).		
Materials indicated on bill of materials present only and labelled appropriately (circle as required): <i>WIP / Approved</i>		
All set-up materials are rejected – in clearly labelled reject containers.		
Ensure recipe version(s) in use for B1240 and Visio Chrom are in compliance with section 1.1 of the PI.		
Sample Blisters: (2 sample blisters – one from each lane):		
- Batch Number and Expiry Date are correct. Record: Batch Number <i>LOT 23045</i> Expiry Date: <i>EXP 04/2026</i>		
VERIFY FORMAT/LAYOUT AGAINST BLISTER SPECIFICATION		
- Pockets are aligned properly at sealing tool-check pockets of blister are not deformed.		
- Lid foil is aligned properly with format foil by checking printing on lid foil are in the correct position.		
- Confirmed that leak test performed and within specification (passed).		
- Retention samples of compliant blisters taken and included with batch record (2 sample blisters – one from each lane).		
I declare that the above checks have been completed and are completed to satisfaction.		
LINE CLEARANCE APPROVED BY/DATE: <i>J 08 Sep 23</i>		

Authorisation signature and date:

J 08 May 23

Reference: SOP-269

Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.

TE-015-03

pharmaxis	FM-301-19	Date Effective: 10-May-23	Page 1 of 3
UHLMANN B1240 FORMAT 10 SETUP CHECKLIST			
Issued By / Date (QA): <i>SV 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

PRODUCT DESCRIPTION – to be completed by Production Supervisor/Packaging Leader

Blister Description:	<i>Branchitol Blister Strip Australia and Europe</i>		
Blister Specification Number:	<i>SP-623-08</i>		
Batch Number to be Embossed on Blister:	<i>LOT 23045</i>		
Expiry Date to be Embossed on Blister: (N/A if no expiry on blister)	<i>EXP 07/2026</i>		
VisioChrom Recipe to be used (circle):	<i>VC_Bronchitol 10s</i>	VC_RedOpaque 10s	VC_Demo 10s
Lid Foil Description:	<i>Branchitol Foil Australia and Europe 2022</i>		
Lid Foil Specification Number:	<i>SP-603</i>		
COMPLETED BY/DATE: <i>Prod. Supervisor/Pack. Leader</i>	<i>SV 05 SEP 23</i>		

	Task Description	Completed✓																																				
1.	Trained Operators must issue and weigh lid & form foils (gross weights only required) & record the balance tag- as detailed in SOP-269.	✓																																				
2.	Login to B1240 control system.	✓																																				
3.	Set B1240 to "Format Change-Over" mode.	✓																																				
4.	Load 'Format 10' format recipe into B1240 control system. Record the Version No. and confirm it to be the same as PI Sect 1.1; Version Number loaded: <i>02</i>	✓																																				
5.	Check condition of following tooling and install. <table border="1"> <thead> <tr> <th>Part Name</th> <th>Part Number</th> <th>Completed✓</th> </tr> </thead> <tbody> <tr> <td>Forming Tool</td> <td>59/1646-1A</td> <td>✓</td> </tr> <tr> <td>Sealing Tool</td> <td>59/1646-1A02</td> <td>✓</td> </tr> <tr> <td>Presealing Plate</td> <td>Pharmaxis drg no 097-030 rev2</td> <td>✓</td> </tr> <tr> <td>Punching Tool</td> <td>59/1646-1A</td> <td>✓</td> </tr> <tr> <td>Punching Tool Guide</td> <td>59/1646-1A</td> <td>✓</td> </tr> <tr> <td>SIMTAP Feed Tubes</td> <td>1A (40 off)</td> <td>✓</td> </tr> <tr> <td>SIMTAP Sorting Plate</td> <td>59/1646-1A</td> <td>✓</td> </tr> <tr> <td>SIMTAP Feed Manifold</td> <td>59/1646-1A</td> <td>✓</td> </tr> <tr> <td>Lid Foil Rollers</td> <td>(Various)</td> <td>✓</td> </tr> <tr> <td>Lid Foil Stretching Station</td> <td>(Various)</td> <td>✓</td> </tr> <tr> <td>Form Foil Bush</td> <td>N/A</td> <td>✓</td> </tr> </tbody> </table>	Part Name	Part Number	Completed✓	Forming Tool	59/1646-1A	✓	Sealing Tool	59/1646-1A02	✓	Presealing Plate	Pharmaxis drg no 097-030 rev2	✓	Punching Tool	59/1646-1A	✓	Punching Tool Guide	59/1646-1A	✓	SIMTAP Feed Tubes	1A (40 off)	✓	SIMTAP Sorting Plate	59/1646-1A	✓	SIMTAP Feed Manifold	59/1646-1A	✓	Lid Foil Rollers	(Various)	✓	Lid Foil Stretching Station	(Various)	✓	Form Foil Bush	N/A	✓	
Part Name	Part Number	Completed✓																																				
Forming Tool	59/1646-1A	✓																																				
Sealing Tool	59/1646-1A02	✓																																				
Presealing Plate	Pharmaxis drg no 097-030 rev2	✓																																				
Punching Tool	59/1646-1A	✓																																				
Punching Tool Guide	59/1646-1A	✓																																				
SIMTAP Feed Tubes	1A (40 off)	✓																																				
SIMTAP Sorting Plate	59/1646-1A	✓																																				
SIMTAP Feed Manifold	59/1646-1A	✓																																				
Lid Foil Rollers	(Various)	✓																																				
Lid Foil Stretching Station	(Various)	✓																																				
Form Foil Bush	N/A	✓																																				
6.	Ensure vibrating cylinders are mounted approximately 20° from vertical.	✓																																				
7.	Set SIMTAP switch to '1'.	✓																																				
8.	Wipe lid foil rollers with a sterile wipe dampened with 70% IPA.	✓																																				
9.	Thoroughly wipe foil stretching station clamping tape with a sterile wipe dampened with 70% IPA to remove burrs and/or build-up.	✓																																				

Authorisation signature and date: <i>SV 08 May 23</i>	Reference: SOP-270
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-301-19	Date Effective: 10-May-23	Page 2 of 3
UHLMANN B1240 FORMAT 10 SETUP CHECKLIST			
Issued By / Date (QA): <i>JS 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

	Task Description	Completed (✓)																																																																																	
10.	Set positions of equipment to the following settings.																																																																																		
	<table border="1"> <thead> <tr> <th>Equipment</th> <th>Machine Setting</th> <th>Actual Setting</th> </tr> </thead> <tbody> <tr> <td>Form Foil Reel Stop</td> <td>65 mm</td> <td><i>65 mm</i></td> </tr> <tr> <td>Lid Foil Reel Stop</td> <td>67 mm</td> <td><i>67 mm</i></td> </tr> <tr> <td>Reversing Station</td> <td>145 mm</td> <td><i>145 mm</i></td> </tr> <tr> <td>Compensation Pendulum</td> <td>10 mm (adjustable)</td> <td><i>10 mm</i></td> </tr> <tr> <td>Print Registration</td> <td>33 mm</td> <td><i>33 mm</i></td> </tr> <tr> <td>Registration Sensor Traverse</td> <td>29 mm</td> <td><i>29 mm</i></td> </tr> <tr> <td>Lid Foil Splice Sensor</td> <td>107 mm (adjustable)</td> <td><i>107 mm</i></td> </tr> <tr> <td>Punching Tool</td> <td>926.5 mm (adjustable)</td> <td><i>926.5 mm</i></td> </tr> <tr> <td>Camera Position</td> <td>40 mm (adjustable)</td> <td><i>40 mm</i></td> </tr> <tr> <td>SIMTAP Sorting Plate Stroke Speed</td> <td>480</td> <td><i>480</i></td> </tr> <tr> <td>SIMTAP Base Position</td> <td>18.3 cm</td> <td><i>18.3 cm</i></td> </tr> <tr> <td>SIMTAP Traverse Position</td> <td>54.5 mm</td> <td><i>54.5 mm</i></td> </tr> <tr> <td>SIMTAP Longitudinal Position</td> <td>40.0±3.0 mm</td> <td><i>40.0 mm</i></td> </tr> <tr> <td>SIMTAP Hopper Outlet</td> <td>15 mm</td> <td><i>15 mm</i></td> </tr> <tr> <td>SIMTAP Level Sensor</td> <td>50 mm (left and right)</td> <td><i>52 mm left 52 mm right</i></td> </tr> <tr> <td>SIMTAP Sorting Plate Vibration Pressure</td> <td>2.4 bar (machine stationary)</td> <td><i>2.4 bar</i></td> </tr> <tr> <td>SIMTAP Filling Tube Vibration Pressure</td> <td>3.2 bar (machine stationary)</td> <td><i>3.2 bar</i></td> </tr> <tr> <td colspan="2">Shaded cells indicate critical parameters.</td> <td></td> </tr> <tr> <td>11.</td> <td>Set embossing type, lot number & expiry date VERIFY FORMAT/LAYOUT AGAINST BLISTER SPECIFICATION</td> <td>✓</td> </tr> <tr> <td>12.</td> <td>Feed form foil and lid foil.</td> <td>✓</td> </tr> <tr> <td>13.</td> <td>Set B1240 to "Test Run" mode.</td> <td>✓</td> </tr> <tr> <td>14.</td> <td>Install SIMTAP feed manifold.</td> <td>✓</td> </tr> <tr> <td>15.</td> <td>Start B1240 and check sealing station aligned with blister pockets. Adjust form foil position and compensation pendulum as necessary. Check depth of seal knurling. Adjust if required – refer SOP-270</td> <td>✓</td> </tr> <tr> <td>16.</td> <td>Position guide before pinhole detector.</td> <td>✓</td> </tr> <tr> <td>17.</td> <td>Login to VisioChrom control system. Load VisioChrom Camera control system recipe as per the product description. Record the Recipe Name & Version No loaded and confirm it to be the same as PI Sect 1.1: VisioChrom Control System Recipe: <i>VC - Drumblat 10's</i> Version Number: <i>09</i></td> <td>✓</td> </tr> <tr> <td>18.</td> <td>Check alignment of seal and pockets by inspecting blisters. Adjust alignment by releasing indexing clamps after sealing station.</td> <td>✓</td> </tr> </tbody> </table>		Equipment	Machine Setting	Actual Setting	Form Foil Reel Stop	65 mm	<i>65 mm</i>	Lid Foil Reel Stop	67 mm	<i>67 mm</i>	Reversing Station	145 mm	<i>145 mm</i>	Compensation Pendulum	10 mm (adjustable)	<i>10 mm</i>	Print Registration	33 mm	<i>33 mm</i>	Registration Sensor Traverse	29 mm	<i>29 mm</i>	Lid Foil Splice Sensor	107 mm (adjustable)	<i>107 mm</i>	Punching Tool	926.5 mm (adjustable)	<i>926.5 mm</i>	Camera Position	40 mm (adjustable)	<i>40 mm</i>	SIMTAP Sorting Plate Stroke Speed	480	<i>480</i>	SIMTAP Base Position	18.3 cm	<i>18.3 cm</i>	SIMTAP Traverse Position	54.5 mm	<i>54.5 mm</i>	SIMTAP Longitudinal Position	40.0±3.0 mm	<i>40.0 mm</i>	SIMTAP Hopper Outlet	15 mm	<i>15 mm</i>	SIMTAP Level Sensor	50 mm (left and right)	<i>52 mm left 52 mm right</i>	SIMTAP Sorting Plate Vibration Pressure	2.4 bar (machine stationary)	<i>2.4 bar</i>	SIMTAP Filling Tube Vibration Pressure	3.2 bar (machine stationary)	<i>3.2 bar</i>	Shaded cells indicate critical parameters.			11.	Set embossing type, lot number & expiry date VERIFY FORMAT/LAYOUT AGAINST BLISTER SPECIFICATION	✓	12.	Feed form foil and lid foil.	✓	13.	Set B1240 to "Test Run" mode.	✓	14.	Install SIMTAP feed manifold.	✓	15.	Start B1240 and check sealing station aligned with blister pockets. Adjust form foil position and compensation pendulum as necessary. Check depth of seal knurling. Adjust if required – refer SOP-270	✓	16.	Position guide before pinhole detector.	✓	17.	Login to VisioChrom control system. Load VisioChrom Camera control system recipe as per the product description. Record the Recipe Name & Version No loaded and confirm it to be the same as PI Sect 1.1: VisioChrom Control System Recipe: <i>VC - Drumblat 10's</i> Version Number: <i>09</i>	✓	18.	Check alignment of seal and pockets by inspecting blisters. Adjust alignment by releasing indexing clamps after sealing station.	✓
	Equipment		Machine Setting	Actual Setting																																																																															
	Form Foil Reel Stop		65 mm	<i>65 mm</i>																																																																															
	Lid Foil Reel Stop		67 mm	<i>67 mm</i>																																																																															
	Reversing Station		145 mm	<i>145 mm</i>																																																																															
	Compensation Pendulum		10 mm (adjustable)	<i>10 mm</i>																																																																															
	Print Registration		33 mm	<i>33 mm</i>																																																																															
	Registration Sensor Traverse		29 mm	<i>29 mm</i>																																																																															
	Lid Foil Splice Sensor		107 mm (adjustable)	<i>107 mm</i>																																																																															
	Punching Tool		926.5 mm (adjustable)	<i>926.5 mm</i>																																																																															
	Camera Position		40 mm (adjustable)	<i>40 mm</i>																																																																															
	SIMTAP Sorting Plate Stroke Speed		480	<i>480</i>																																																																															
	SIMTAP Base Position		18.3 cm	<i>18.3 cm</i>																																																																															
	SIMTAP Traverse Position		54.5 mm	<i>54.5 mm</i>																																																																															
	SIMTAP Longitudinal Position		40.0±3.0 mm	<i>40.0 mm</i>																																																																															
	SIMTAP Hopper Outlet		15 mm	<i>15 mm</i>																																																																															
	SIMTAP Level Sensor		50 mm (left and right)	<i>52 mm left 52 mm right</i>																																																																															
SIMTAP Sorting Plate Vibration Pressure	2.4 bar (machine stationary)	<i>2.4 bar</i>																																																																																	
SIMTAP Filling Tube Vibration Pressure	3.2 bar (machine stationary)	<i>3.2 bar</i>																																																																																	
Shaded cells indicate critical parameters.																																																																																			
11.	Set embossing type, lot number & expiry date VERIFY FORMAT/LAYOUT AGAINST BLISTER SPECIFICATION	✓																																																																																	
12.	Feed form foil and lid foil.	✓																																																																																	
13.	Set B1240 to "Test Run" mode.	✓																																																																																	
14.	Install SIMTAP feed manifold.	✓																																																																																	
15.	Start B1240 and check sealing station aligned with blister pockets. Adjust form foil position and compensation pendulum as necessary. Check depth of seal knurling. Adjust if required – refer SOP-270	✓																																																																																	
16.	Position guide before pinhole detector.	✓																																																																																	
17.	Login to VisioChrom control system. Load VisioChrom Camera control system recipe as per the product description. Record the Recipe Name & Version No loaded and confirm it to be the same as PI Sect 1.1: VisioChrom Control System Recipe: <i>VC - Drumblat 10's</i> Version Number: <i>09</i>	✓																																																																																	
18.	Check alignment of seal and pockets by inspecting blisters. Adjust alignment by releasing indexing clamps after sealing station.	✓																																																																																	

Authorisation signature and date: <i>[Signature] 08 May 23</i>	Reference: SOP-270
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-301-19	Date Effective: 10-May-23	Page 3 of 3
UHLMANN B1240 FORMAT 10 SETUP CHECKLIST			
Issued By / Date (QA): <i>JK 01 Sep 23</i>		Lot Number: <i>M23-045</i>	

	Task Description	Completed (✓)
19.	Run machine & check lid foil print registration mark at the blister seam. If the registration mark is immediately left & to the side of the seam, the machine will automatically align. Otherwise see SOP-269 for foil adjustment.	✓
20.	Ensure that control block 2/3 pneumatic regulator is set to 4 bar.	✓
21.	Check VisioChrom camera image is aligned with pockets and adjust camera position as necessary.	✓
22.	Check chiller water level is good, pump operating and controlling to setpoint temperature.	✓
23.	Set B1240 to "Production" mode to activate the detection system.	✓
24.	Run B1240 for approximately 100 cycles and check blisters for: <ul style="list-style-type: none"> No visual defects Sealing aligned to pockets Embossing is correct and clear Lid foil printing aligned Sufficient sealing border around perimeter 	✓
25.	Sample preliminary leak test blisters as detailed in SOP-270. Manually sample 6 blisters (3 from each lane) & stop the machine. Send blister samples to QC for seal integrity testing.	✓
26.	Remove two blisters for retention. Label as R1 and R2 appropriately (where 1 and 2 correspond to lane 1 and 2).	✓
27.	Perform pinhole challenge test as per SOP-270.	✓
28.	Raise compensation pendulum (to prevent foil slippage on loss of compressed air).	✓
29.	Logout of B1240 control system and VisioChrom camera if not starting batch immediately after setup.	✓ / NA
30.	Additional Comments: <i>N/A</i>	
31.	Assembly and set-up completed according to SOP-270 by: Print Name: <i>German G.</i> Sign/Date: <i>gg 07 Sep 23</i>	
32.	Blister Seal Integrity Testing Integrity of sealing tested as per TM-046. If blisters fail leak test, refer to SOP-270 to rectify leak. Repeat preliminary leak test sampling. Once a 'PASS' result is achieved, continue batch processing as per PI. Print Name: <i>German G.</i> Sign/Date: <i>gg 08 Sep 23</i> Operator/Maintenance Engineer	

Authorisation signature and date: <i>ASD 08 May 23</i>	Reference: SOP-270
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmoxis	FM-297-09	Date Effective: 28-Aug-20	Page 1 of 2
B1240 UHLMANN BLISTER THERMOFORMER (00407) MAJOR CLEAN CHECKLIST			

Machine Cleaned After Lot:

M23-034

Batch of 70% IPA used:

IPA-0496

Completed By / Date:

gg 06 Sep 23

REFER TO SOP-271 FOR FULL CLEANING PROCEDURE

MACHINE PARTS	
1.	Wipe with dry non-linting wipe only
	Form film splice detection sensor
	Lid Foil splice detection sensor
	VisioChrom camera lens
	COMPLETED BY/ DATE: gg 06 Sep 23
2.	Wipe with non-linting wipe dampened with 70% IPA
	Forming Tools
	Punching Tool
	COMPLETED BY/ DATE: gg 06 Sep 23
3.	Flush and rinse with WFI; Spray with 70% IPA
	Shut-Off Housing Dosing Rings
	Chip Collectors
	Vibratory Chute
	Sorting Plate and Frame
	Fill Tubes - gently clean internal surfaces with bottle brush.
	Product Chute Covers
	SIMTAP Manifold (disassemble)
	Feeding Hoses
	Perforated Plates
	Hoppers
	COMPLETED BY/ DATE: gg 06 Sep 23
4.	Clean with brass brush; Clear out the film residue using a non-linting wipe; Wipe with non-linting wipe dampened with 70% IPA
	Sealing Plates & Product Haloes
	COMPLETED BY/ DATE: gg 06 Sep 23
5.	Clean with WFI/ Wipe with non-linting wipe dampened with 70% IPA
	Control Panel
	Machine Base
	Blister Slide
	Punching Tool Conveyor
	Form Foil Bush
	(Note: Check for any wear/tear on the bush, any issues should be notified to production supervisor immediately.)
	Splice Table
COMPLETED BY/ DATE: gg 06 Sep 23	
6.	Move the cleaned & dried sanitised parts to Clean Equipment Store Room (B1.026).
	COMPLETED BY/ DATE: gg 06 Sep 23

Authorisation signature and date:


gg 20 Aug 20


Reference: SOP-271

Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.

TE-002-03

pharmaxis	FM-297-09	Date Effective: 28-Aug-20	Page 2 of 2
B1240 UHLMANN BLISTER THERMOFORMER (00407) MAJOR CLEAN CHECKLIST			

EQUIPMENT INSPECTION & APPROVAL	
7.	A second person (trained in TD-038: Responsibility of a Checker) must inspect and approve (if satisfactory) the cleaning of B1240 Blister machine (tag no. 00407) and all the parts. The Checker is to inspect / check the following:
	The equipment is visually clean, dry and fit for use.
	Equipment parts are visually clean and dry.
	Thoroughly check equipment for visual signs of any foreign particles (i.e. dust, rust on any part, etc.) or any visible equipment or gasket damage.
	External surfaces of the equipment are clean and dry.
	Area surrounding the equipment is clean, dry and free of rubbish.
	Ensure that cleaning tags are attached to the major equipment and its cleaned detachable parts as per SOP-191.
	Sign below only if the checker is fully satisfied that the equipment is clean and fit for use. Report any non-conformities to the Production Supervisor:
Checked By/ Date: <u> 06 Sep 23</u>	

Authorisation signature and date: <u> 20 Aug 20</u>	Reference: SOP-271
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-231-04	Date Effective: 25-Jul-18	Page 1 of 1
CLEANING TAG/EQUIPMENT STATUS			

EQUIPMENT

Equipment Name	Conveyor Belt
Equipment Tag Number	00910

Cleaning Status:

Type of clean	Completed By	Date
Minor Clean	<i>~</i>	07 Sep 23
Minor Clean Awaiting Major Clean		
Major Clean		N/A 07 Sep 23
Major Clean Approved		<i>~</i>
Resanitisation? (Include sign and date)		


Authorisation signature and date: <i>Steph Dally 12 Feb 18</i>	Reference: SOP-142
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	


pharmaxis	FM-231-04	Date Effective: 25-Jul-18	Page 1 of 1
CLEANING TAG/EQUIPMENT STATUS			

EQUIPMENT

Equipment Name	WILMANN BIZ4D (Bronchitol 10's)
Equipment Tag Number	00407

Cleaning Status:

Type of clean	Completed By	Date
Minor Clean	N/A	N/A
Minor Clean Awaiting Major Clean	N/A	N/A
Major Clean	gg	06 Sep 23
Major Clean Approved		06 Sep 23
Resanitisation? (Include sign and date)	N/A	N/A

Authorisation signature and date:  12 Jul 18	Reference: SOP-142
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-231-04	Date Effective: 25-Jul-18	Page 1 of 1
CLEANING TAG/EQUIPMENT STATUS			

EQUIPMENT

Equipment Name	WILMANN B1240 (Parts)
Equipment Tag Number	00407

Cleaning Status:

Type of clean	Completed By	Date
Minor Clean	N/A	N/A
Minor Clean Awaiting Major Clean	N/A	N/A
Major Clean	gg	06 SEP 23
Major Clean Approved	[Signature]	06 Sep 23
Resanitisation? (Include sign and date)	N/A	* N/A

* Amended on DAY 105 OCT 23

Authorisation signature and date: [Signature] 12 Feb 18	Reference: SOP-142
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10

10-1-10


10-1-10


pharmaxis	FM-231-04	Date Effective: 25-Jul-18	Page 1 of 1
CLEANING TAG/EQUIPMENT STATUS			

EQUIPMENT

Equipment Name	UHLMANN B124D
Equipment Tag Number	00407

Cleaning Status:

Type of clean	Completed By	Date
Minor Clean	N/A	N/A
Minor Clean Awaiting Major Clean	N/A	N/A
Major Clean	gg	06 Sep 23
Major Clean Approved		06 Sep 23
Resanitisation? (Include sign and date)	N/A	N/A

Authorisation signature and date:  12 Feb 18	Reference: SOP-142
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

1944

pharmaxis	FM-305-09	Date Effective: 13-Jun-23	Page 1 of 4
OUTPUT BLISTERS FROM UHLMANN B1240			
Lot Number: M23-045		Issued By/Date: <i>JK</i> 01 Sep 23	

SECTION 1 - to be completed by Production Supervisor/Packaging Leader/Delegate

Lot Number	M23-045	Expiry Date	EXP 07/2026
Blister Description	Branchitel Bister Strip Au EU	Specification Number (include version)	SP-623-08
Completed By/Date: <i>[Signature]</i> 1 05 SEP 23			

SECTION 2 – Perform room clearance.

Room clearance must be performed by an appropriately trained person.	✓ each step once complete
Correct area (Secondary Packaging – B1.034) is in use.	✓
The outer door of the Secondary Packaging room is status labelled.	✓
All the documentation, products and materials from the previous production operation have been removed from the Secondary Packaging room.	✓
Ensure that the appropriate cleaning has been performed/approved for the Conveyor Belt (00910).	✓
Cleaning tag (FM-231) for Conveyor Belt (00910) has been attached to this FM-305.	✓
LB-180 (Secondary Packaging) and LB-207 (Conveyor Belt) are checked for completeness (eg., Login, batch details, minor clean/ daily clean performed).	✓
No alarms on EMS system (ensure alarm and strobe lighting not active).	✓
Comments: B1.034	
I declare that the above checks have been completed to my satisfaction.	
ROOM CLEARANCE APPROVED BY/DATE: <i>JK</i> 05 SEP 23	

Sign & Date in the Secondary Packaging Area (LB-180) logbook upon successful completion of the Room Clearance.

Authorisation signature and date: <i>JK</i> 13 JUN 23	Reference: SOP-297
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-305-09	Date Effective: 13-Jun-23	Page 2 of 4
OUTPUT BLISTERS FROM UHLMANN B1240			
Lot Number: M23-045		Issued By/Date: <i>[Signature]</i> 01 Sep 23	

SECTION 3- PROCEDURE

1. Use shippers 52 x 36 x 47cm for collecting blisters.
2. Pack blisters in configuration described on page 4 of this form.
3. Collect the blisters from conveyor belt, count and pack the blisters in to the shipper as detailed in SOP-297. Separate each layer of blister with blank paper sheet (A3)/inserts.
4. Once shipper is full, complete WIP label (CL-008). The shippers are numbered sequentially. Complete the table below after each shipper is filled/complete.

Shipper No.	Date	Shipper End Time	Blister Quantity Packed per shipper	Initial
1	08 Sep 23	12:58	1700	<i>[Initial]</i>
2	08 Sep 23	13:24	1700	<i>[Initial]</i>
3	11 Sep 23	07:08	1700	<i>[Initial]</i>
4	11 Sep 23	07:36	1700	<i>[Initial]</i>
5	11 Sep 23	07:58	1700	<i>[Initial]</i>
6	11 Sep 23	08:22	1700	<i>[Initial]</i>
7	11 Sep 23	09:35	1700	<i>[Initial]</i>
8	11 Sep 23	10:00	1700	<i>[Initial]</i>
9	11 Sep 23	10:35	1700	<i>[Initial]</i>
10	11 Sep 23	11:05	1700	<i>[Initial]</i>
11	11 Sep 23	12:13	1700	<i>[Initial]</i>
12	11 Sep 23	12:39	1700	<i>[Initial]</i>
13	11 Sep 23	13:13	1700	<i>[Initial]</i>
14	11 Sep 23	13:36	1700	<i>[Initial]</i>
15	12 Sep 23	06:59	1700	<i>[Initial]</i>
16	12 Sep 23	07:35	1700	<i>[Initial]</i>
17	12 Sep 23	07:51	656	<i>[Initial]</i>
18				
19				
20				
21				
22				
23				
24				

Authorisation signature and date: <i>[Signature]</i> 13 JUN 23	Reference: SOP-297
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

pharmaxis	FM-305-09	Date Effective: 13-Jun-23	Page 3 of 4
OUTPUT BLISTERS FROM UHLMANN B1240			
Lot Number: M23-045		Issued By/Date: <i>[Signature]</i> 01 Sep 23	

Shipper No.	Date	Shipper End Time	Blister Quantity Packed per shipper	Initial
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				


SECTION 3. RECONCILIATION

Total number of blisters rejected during the collection: <u>1</u>
Completed By/Date: <u><i>[Signature]</i> 12 Sep 23</u>
Ensure that the data recorded on all the shipper labels is accurate.
Checked By/Date: <u><i>[Signature]</i> 12 SEP 23</u>
PRODUCTION DOCUMENT REVIEW: To be completed by production delegate
Checked By/Date: <u><i>[Signature]</i> 18 SEP 23</u>

Authorisation signature and date: <i>[Signature]</i> 13 JUN 23	Reference: SOP-297
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	


pharmaxis	FM-305-09	Date Effective: 13-Jun-23	Page 4 of 4
OUTPUT BLISTERS FROM UHLMANN B1240			
Lot Number: M23-045		Issued By/Date: <i>[Signature]</i> 01 Sep 23	

Aridol /Osmohale Blister Strip 1




130 Blisters per row
 5 rows = 5 x 130 blisters = 650 blisters
 4 Layers = 4 x 650 blisters = **2600 blisters**

Aridol/Osmohale Blister Strip 2 and 3



85 Blisters per row
 5 rows = 5 x 85 blisters = 425 blisters
 4 Layers = 4 x 425 blisters = **1700 blisters**

Bronchitol Blisters (applies to any Format 10 blister)



85 Blisters per row
 5 rows = 5 x 85 blisters = 425 blisters
 4 Layers = 4 x 425 blisters = **1700 blisters**

Authorisation signature and date: <i>[Signature]</i> 13 JUN 23	Reference: SOP-297
Confidential. Staff must comply with the current version of the document. Reproductions of this document are uncontrolled and may not be the current version. Not for unauthorised distribution.	
TE-002-03	

Statistic - Thermoformer



Uhlmann

Format 4 Format 10-02

M23-045 Bronchitol Blister Strip An/Eu

Good blister	28134
Misfilled blisters	4120
MSR init	232
Total rejection	843
Missing product	855
Empty blisters	1264
Test run	0
Splice forming film	42
Overheated web heating station	0
Forming fault cause emergency stop	28
Pinhole detection	12
Fill control not ready	0
Error Data-Valid fill control	0
Splice lid foil	6
Code error	0
Overheated web sealing station	838
Overload sealing station	0

gg 12 Sep 23


 UNIVERSITÄT
 AN DER SAAR

Fachbereich Wirtschaftswissenschaften
 Fachhochschule für Technik und Wirtschaft

Fachbereich Wirtschaftswissenschaften
 Fachhochschule für Technik und Wirtschaft

1. Semester	1. Semester
2. Semester	2. Semester
3. Semester	3. Semester
4. Semester	4. Semester
5. Semester	5. Semester
6. Semester	6. Semester
7. Semester	7. Semester
8. Semester	8. Semester
9. Semester	9. Semester
10. Semester	10. Semester
11. Semester	11. Semester
12. Semester	12. Semester
13. Semester	13. Semester
14. Semester	14. Semester
15. Semester	15. Semester
16. Semester	16. Semester
17. Semester	17. Semester
18. Semester	18. Semester
19. Semester	19. Semester
20. Semester	20. Semester
21. Semester	21. Semester
22. Semester	22. Semester
23. Semester	23. Semester
24. Semester	24. Semester
25. Semester	25. Semester
26. Semester	26. Semester
27. Semester	27. Semester
28. Semester	28. Semester
29. Semester	29. Semester
30. Semester	30. Semester
31. Semester	31. Semester
32. Semester	32. Semester
33. Semester	33. Semester
34. Semester	34. Semester
35. Semester	35. Semester
36. Semester	36. Semester
37. Semester	37. Semester
38. Semester	38. Semester
39. Semester	39. Semester
40. Semester	40. Semester
41. Semester	41. Semester
42. Semester	42. Semester
43. Semester	43. Semester
44. Semester	44. Semester
45. Semester	45. Semester
46. Semester	46. Semester
47. Semester	47. Semester
48. Semester	48. Semester
49. Semester	49. Semester
50. Semester	50. Semester
51. Semester	51. Semester
52. Semester	52. Semester
53. Semester	53. Semester
54. Semester	54. Semester
55. Semester	55. Semester
56. Semester	56. Semester
57. Semester	57. Semester
58. Semester	58. Semester
59. Semester	59. Semester
60. Semester	60. Semester
61. Semester	61. Semester
62. Semester	62. Semester
63. Semester	63. Semester
64. Semester	64. Semester
65. Semester	65. Semester
66. Semester	66. Semester
67. Semester	67. Semester
68. Semester	68. Semester
69. Semester	69. Semester
70. Semester	70. Semester
71. Semester	71. Semester
72. Semester	72. Semester
73. Semester	73. Semester
74. Semester	74. Semester
75. Semester	75. Semester
76. Semester	76. Semester
77. Semester	77. Semester
78. Semester	78. Semester
79. Semester	79. Semester
80. Semester	80. Semester
81. Semester	81. Semester
82. Semester	82. Semester
83. Semester	83. Semester
84. Semester	84. Semester
85. Semester	85. Semester
86. Semester	86. Semester
87. Semester	87. Semester
88. Semester	88. Semester
89. Semester	89. Semester
90. Semester	90. Semester
91. Semester	91. Semester
92. Semester	92. Semester
93. Semester	93. Semester
94. Semester	94. Semester
95. Semester	95. Semester
96. Semester	96. Semester
97. Semester	97. Semester
98. Semester	98. Semester
99. Semester	99. Semester
100. Semester	100. Semester

Fachbereich Wirtschaftswissenschaften
 Fachhochschule für Technik und Wirtschaft

Fachbereich Wirtschaftswissenschaften
 Fachhochschule für Technik und Wirtschaft



31. August 2023

Production Order - Job Card

Pharmaxis Master Data

Page 1

RowenaG

Production Order

Prod. Order No. 23/109
Status Released

Item Bronchitol Blister AU EU
Item No. 200045
SP No. SP-623-08
Lot No. M23-045

Starting Date 24/12/18
Starting Time 12:00:00 PM
Ending Date 31/12/18
Ending Time 2:00:00 PM

Expected Quantity 40,000 Each
Output Bin Code 20-F1K PACK. STAGGIN

Routing

Routing No. SP-623

Operation No. 10
Type Work Center
No. B1.022
Description Primary Packaging Room #20
Time Needed 42 HR

Material Requirements

Item No.	SP No.	Item	Expected Quantity	Unit of Measure Code	Consumed Quantity	Lot No.	Bin Code
200042	SP-501-10	40mg Print Clear capsule #20RR	36,452	G			
300144	SP-603-06	Bronchitol Foil AU/EU - #20RR	29	KG			
300153	SP-627-04	Blister Form Foil 160mm -20RR	96	KG			

Output

Gross Output _____ Each
Samples _____ Each
Net Output _____ Each